



artacent[®]

flex

High-quality, regenerative amniotic tissue graft that is dual-layered, creating an ideal gliding surface for adjacent tissues.

to us, it's personal.

tides medical[®] stands behind our advanced products and expert service, so you can...

Help patients get the treatments they deserve. By providing high-quality, safe products for patients and supporting clinicians throughout the process, we ensure more people can benefit from advanced treatment.

INCREDIBLY VERSATILE PRODUCT THAT CAN BE USED IN VARIOUS SURGICAL PROCEDURES:

- Orthopaedic
- Spinal
- Urological
- Vascular
- Ophthalmic
- Plastics
- Colorectal
- General Surgery

888-494-4441

www.tidesmedical.com



tidesmedical[®]

advanced products. expert service.

artacent[®]

flex

High-quality, regenerative amniotic tissue graft that is dual-layered, creating an ideal gliding surface for adjacent tissues.



Artacent Flex[®] is a safe, natural biologic barrier that maintains the original powerful characteristics of amniotic tissue. The unique design of Artacent Flex allows for easy manipulation and repositioning, making it a flexible, dependable option for a variety of surgical applications.

EASE-OF-USE: Can be placed on either side facing the wound, remains in place and doesn't require sutures.


FAVORABLE: Ample inventory options and overnight shipping, favorable customer programs, with minimal risk

NATURAL: 100 % pure amniotic tissue processed using Tides' proprietary Artacleanse[®] process.

FLEXIBLE: Stores dry at room temperature with no need for refrigeration. Product remains safe and effective with a shelf life of 5 years.

ARTACENT FLEX PART NUMBERS AND SIZES

Part Number	Size
1. ART1101	1 x 1 cm
2. ART1202	2 x 2 cm
3. ART1203	2 x 3 cm
4. ART1404	4 x 4 cm
5. ART1406	4 x 6 cm
6. ART1408	4 x 8 cm
7. ART1016	16 mm disk



Contact the sales and service professionals in our nationwide network to learn how we take your business personally.

tidesmedical

advanced products. expert service.

888-494-4441

www.tidesmedical.com